


## Variability of Affinity Selection



Standard Curve of Hu Transferrin						
Amount (ng)	Run 1	Run 2	Run 3	Average	Stdev	CV
5000	210.1	225.8	222.8	219.6	8.3	3.8
2500	202.6	211.0	218.8	210.8	8.1	3.8
1250	175.9	178.1	179.5	177.8	1.8	1.0
500	131.6	130.5	130.9	131.0	0.6	0.4
75	62.7	62.6	64.4	63.2	1.0	1.6
10	28.7	30.2	27.9	28.9	1.2	4.0